

## d-PETT<sup>ultra</sup> Single and Multi-Channel Pipette



- Magnet assisted piston = high precision + reduced pipetting force
- Ergonomic design for easy pipetting operation
- Fatigue free working with convenient finger rest
- Fully autoclavable & UV resistance
- 4-digit volume counter for a more volume setting
- Tip ejector with built-in shock absorber for reduced thumb impact
- Corrosion free stainless steel & fully autoclavable tip ejector
- Color coding for easy identification
- Easy cleaning & in-house calibration
- Calibrated in ISO 17025 (NABL) accredited laboratory. Supplied with individual QC & calibration report according to ISO 8655-5
- Min Vol - 0.1 µl & Max Vol - 10,000 µl
- 3-Year warranty

Buy Any 5x d-PETT<sup>ultra</sup> and  
get one Beautiful d-  
PETT<sup>ultra</sup> Pipette Carousel  
Stand Free



Variable Volume Micropipette		Fixed Volume Micropipette	
Cat. No.	Volume Range [ $\mu$ l ]	Cat. No.	Volume Range [ $\mu$ l ]
DPU-2.5	0.1-2.5	DPU-205	2.5
DPU-10	0.5-10	DPU-5	5
DPU-20	2-20	DPU-10	10
DPU-50	5-50	DPU-20	20
DPU-100	10-100	DPU-25	25
		DPU-50	50
DPU-200	20-200	DPU-100	100
		DPU-200	200
DPU-1000	100-1000	DPU-250	250
		DPU-500	500
DPU-5000	500-5000	DPU-1000	1000
		DPU-2000	2000
DPU-10000	1000-10000	DPU-5000	5000
		DPU-10000	10000

8 Channel Pipette		12 Channel Pipette	
Cat. No.	Volume Range [ $\mu$ l ]	Cat. No.	Volume Range [ $\mu$ l ]
DPU8M-10	0.5 - 10 $\mu$ l	DPU12M-10	0.5 - 10 $\mu$ l
DPU8M -20	2 - 20 $\mu$ l	DPU12M-20	2 - 20 $\mu$ l
DPU8M -50	5- 50 $\mu$ l	DPU12M -50	5 - 50 $\mu$ l
DPU8M -100	10 - 100 $\mu$ l	DPU12M -100	10 - 100 $\mu$ l
DPU8M -200	20 - 200 $\mu$ l	DPU12M -200	20 - 200 $\mu$ l
DPU8M -300	30 - 300 $\mu$ l	DPU12M -300	30 - 300 $\mu$ l



# 3 YEAR WARRANTY

All pipettes carry three years warranty from the date of purchase, against defects in material and workmanship. We will repair or replace the defective product upon prompt notification in compliance with the conditions listed below

The warranty against defects in workmanship or materials will be accepted as long as the product has been used in compliance with the instructions for care and use. The warranty will be considered void if pipette is exposed to physical or chemical abuse. Routine cleaning and recalibration are not covered under the terms of the warranty.

This warranty is subject to the following conditions:

- The customer must notify its authorized distributor in writing immediately after sighting of the defect.
- The product must be returned to its authorized distributor upon request. The complete cost of forwarding in such cases including transportation, applicable Government taxes, insurance is to be borne by the customer.
- Our responsibility extends only to product defects. The warranty does not apply if, in opinion of the product has been damaged by accident, misuse, chemical or physical corrosion and/or by service or modification by someone other than an authorized service provider.

# CERTIFICATE OF CONFORMITY

We declare that all pipettes are manufactured and verified by technicians authorized by us under strictly defined and monitored conditions, according to the ISO 8655-6 standard.

Before dispatching, all pipettes are calibrated in the laboratory by using distilled water grade 3 (ISO 3696) at controlled room temperature, relative air humidity (minimum 50%) and barometric pressure (1013hPa).

Ten measurements are taken at each test volume (nominal, 50% nominal, minimum or 10% nominal whichever is higher). Balances are regularly checked with weight traceable to national and international standards.

Quality system and procedure is well defined and maintained as per ISO 9001:2015 and is regularly audited by external agencies.